

## Haematology

Part no.	Description	unit
----------	-------------	------

### Instruments

HI5110	<b>AMP Accos 5110 - automated 5-part differential hematology analyzer 27 parameters, 2 histo-, 2 scatter-, 2 3D-stereograms, 110 samples / hour, incl. rack sampler, PC, monitor, keyboard, mouse, printer, software, accessories and operation manual</b>	1 pc.
	<i>Accessories included:</i>	
	Sample Tube Rack	12 pcs.
	Reagent Tubing, color-coded	4 pcs.
	Waste Tubing	1 pc.
	Waste Detection Sensor	1 pc.
	Waste Container	2 pcs.
	Peristaltic Pump Tubing	1 pc.
	Filter	2 pcs.
	Diluent Syringe Sealing	2 pcs.
	Lyse Syringe Sealing	2 pcs.
	Grease	1 pc.
	Probe Cleaning Solution	1 pc.
	Serial Port Interface Cable	1 pc.
	Fuse (3.15A / 250V)	2 pcs.
	Power Cord (european standard)	4 pcs.
	Grounding Wire	1 pc.

### Product Documentation

HM5110	Operation Manual - AMP Accos 5110, English	1 pc.
HM9501	Service Manual - AMP Accos 5110	1 pc.

### Reagents

HR5020	AMP HemoDil AS 5	20 ltr.
HR5101	AMP HemoLyse AS 5	1 ltr.
HR5110	AMP HemoSheath AS	10 ltr.
HR5120	AMP HemoTerge AS 5	20 ltr.
HR5205	AMP HemoStain R	50 x 3.7 mL

### Control Blood

HR9352	AMP HemoTrol 5D - Kit S	1 H, 2 N, 1 L x 3 mL
HR9353	AMP HemoTrol 5D - Kit M	2 H, 4 N, 2 L x 3 mL
HR9354	AMP HemoTrol 5D - Kit L	3 H, 6 N, 3 L x 3 mL